WARNING NOTICE: The experiments described in these materials are potentially hazardous and require a high level of safety training, special facilities and equipment, and supervision by appropriate individuals. You bear the sole responsibility, liability, and risk for the implementation of such safety procedures and measures. MIT shall have no responsibility, liability, or risk for the content or implementation of any of the material presented. Legal Notices

BRADFORD ASSAY FOR DETERMINING PROTEIN CONCENTRATION

All preparations should be done in/with plastic because the Bradford Reagent reacts with glass.

- 1. Dilute protein assay reagent (Biorad) 1:4 with milli-Q water (usually 5 mL of reagent and 15 mL of water is enough).
- 2. Add 1 mL of reagent to 200 uL of each protein sample.
- 3. Measure the absorbance of each sample at 595 nm.

To determine the concentration of your sample, plot a standard curve and extrapolate back to your reading. Bovine serum albumin (BSA) is usually the protein standard of choice. The table below indicates the amounts of 400 ug/mL BSA you should mix with water or buffer to give the concentrations shown.

Concentration	Amount of BSA (uL)	Amount of Water/Buffer (uL)
250	125	75
200	100	100
150	75	125
100	50	150
50	25	175
25	12.5	187.5
0	0	200