WARNING NOTICE: The experiments described in these materials are potentially hazardous and require a high level of safety training, special facilities and equipment, and supervision by appropriate individuals. You bear the sole responsibility, liability, and risk for the implementation of such safety procedures and measures. MIT shall have no responsibility, liability, or risk for the content or implementation of any of the material presented. Legal Notices

## ETHANOL PRECIPITATION OF NUCLEIC ACIDS IN AQUEOUS SOLUTION

- 1. Determine approximate volume of solution
- Add either 1/2 volume of 7.5 M ammonium acetate -or-1/10 volume of 3 M sodium acetate pH 5.3
- 3. Add 2.5 to 3 (original) volumes of cold ethanol
- 4. Mix, and incubate on ice (or at -20°C) 15 minutes to overnight
- 5. Centrifuge 15-20 minutes at high speed (in microfuge for small volumes, or in SS34 at 8 000 rpm)
- 6. Discard supernatant
- 7. Wash pellet by adding one original volume cold 70% ethanol
- 8. Recentrifuge if pellet has become dislodged
- 9. Discard supernantant
- 10. Dry pellet either by air drying at room temperature or with Speed-Vac